

Lab ARMOR[®] DRYTEMP[™]

The DryTemp, Bead Boxes, and Beads combine for space efficiency.



The DryTemp digital dry bath offers control over both temperature and time, eliminating the need for an external thermometer or timer. Using the Lab Armor Bead Boxes filled with Beads in the DryTemp instead of the aluminum drilled blocks allows for testing of any size tube. The DryTemp has a built-in user recalibration function, to guarantee long term temperature accuracy.



- ✓ **Control:** Digital accuracy and convenience.
- ✓ **Removable Bead Boxes:** Keep samples heated away from DryTemp unit.
- ✓ **Stays Clean:** Bead Box stays dry and protected from water-borne microbial contaminants.
- ✓ **Organized:** Beads hold capped and non-watertight vessels in place without racks or accessories.
- ✓ **Energy Efficient:** Beads have low biocide use and are recyclable.

Lab Armor[®] DryTemp[™] Baths are durable, dependable, and are available with a variety of safety and convenience features:

- ✓ Digital microprocessor control.
- ✓ Compatible with 2 Single or 1 Double Bead Box.
- ✓ Timer function.
- ✓ 8.7 x 10.3 x 3.2 in (22.1 x 26.2 x 8.1 cm).
- ✓ Temperature Range: 5°C above ambient to 150°C.
- ✓ Electrical Specifications: Volts: 120V or 230, Hertz: 50/60Hz Watts: 500W, Amps: 5.0A.
- ✓ 2-year warranty.

Get a DryTemp[™] and upgrade your entire lab experience.

Part ID	Description
DryTemp™	
52120-200	DryTemp Digital Dry Bath 120 V
52230-200	DryTemp Digital Dry Bath 230 V
Bead Box™	
5620753	Single Box Empty
5620752	Double Box Empty
52300-BLK	Single Box w/ 0.25 liters of Lab Armor Beads
52400-BLK	Double Box w/ 0.50 liters of Lab Armor Beads
DryTemp™ + Bead Box™	
52120-205	120 Volt DryTemp w/ Double Box and 0.50 liters of Lab Armor Beads
52120-210	120 Volt DryTemp w/ two Single Boxes and 0.50 liters of Lab Armor Beads
52230-205	230 Volt DryTemp w/ Double Box and 0.50 liters of Lab Armor Beads
52230-210	230 Volt DryTemp w/ two Single Boxes and 0.50 liters of Lab Armor Beads



Technical Support

For additional product and technical information, such as product manuals or technical articles on the use of DryTemp, Bead Boxes, and Lab Armor Beads for common laboratory applications, please visit labarmor.com. For further assistance, please email our Technical Support team at sales@labarmor.com.

www.labarmor.com | sales@labarmor.com | 1-469-496-2350

Notifications

The trademarks mentioned herein are the property of Lab Armor, LLC.
 Lab Armor® Beads are patented. For laboratory use only. Not for personal use.
 © 2010-2024 Lab Armor. All rights reserved.